

# Surgimpex & Co.

Highest Quality Electrosurgical Instruments





McPherson (Str.)  
10-001-05



McPherson (crv.)  
10-002-05




Jeweler # 5  
10-003-05



Adson  
10-006-10



Jeweler # 7  
10-008-05



Jeweler # 4  
10-007-05

Cat #	Discription	Size
10-001-05	L= 9.0 cm tip	0.5 mm
10-002-05	L= 11.5 cm tip	0.5 mm
10-003-05	L= 11.5 cm tip	0.0 mm
10-006-05	L= 12.1 cm tip	1.0 mm
10-008-05	L= 11.5 cm tip	0.5 mm
10-007-05	L= 11.5 cm tip	0.5 mm





Coaptation  
10-009-05




Coaptation Nadler  
10-010-05



Iris (Str.)  
10-011-10




Iris (crv.)  
10-012-10



Straight  
10-016-10

Straight  
10-018-10

Cat #	Discription	Size
10-009-05	L=10.8 cm	0.5 mm
10-010-05	L=10.2 cm tip	0.5 mm
10-011-10	L=10.8 cm tip	1.0 mm
10-012-10	L=10.8 cm tip	1.0 mm
10-016-10	L=16.0 cm tip	1.0 mm
10-018-10	L=18.0 cm tip	1.0 mm



Cat #	Discription	Size
10-020-10	L=20.0 cm tip	1.0 mm
10-022-10	L=22.0 cm tip	1.0 mm
10-026-10	L=16.0 cm tip	1.0 mm
10-028-10	L=18.0 cm tip	1.0 mm
10-030-10	L=20.0 cm tip	1.0 mm
10-032-10	L=22.0 cm tip	1.0 mm
10-036-10	L=16.0 cm tip	1.0 mm
10-038-10	L=18.0 cm tip	1.0 mm
10-040-10	L=20.0 cm tip	1.0 mm



## European Bipolar Forceps



Bayonet (Ang.Up)  
10-046-10

Bayonet (Ang.Up)  
10-048-10

Bayonet (Ang.Up)  
10-050-10

Bayonet (Ang.Down)  
10-056-10

Bayonet (Ang.Up)  
10-058-10


Bayonet (Ang.Up)  
10-060-10

Cat #	Discription	Size
10-046-10	L=16.0 cm tip	1.0 mm
10-048-10	L=18.0 cm tip	1.0 mm
10-050-10	L=20.0 cm tip	1.0 mm
10-056-10	L=16.0 cm tip	1.0 mm
10-058-10	L=18.0 cm tip	1.0 mm
10-060-10	L=20.0 cm tip	1.0 mm





## European Bipolar Forceps



McPherson (Str.)  
101-001-05

McPherson (crv.)  
101-002-05

Jeweler # 5  
101-003-05

Jeweler # 5  
101-003-05

Jeweler # 7  
101-008-05

Adson  
101-006-10



Straight  
101-016-10



Straight  
101-018-10



Straight  
101-020-10

Cat #	Description	Size
101-001-05	L= 9.0 cm tip	0.5 mm
101-002-05	L= 11.5 cm tip	0.5 mm
101-003-05	L= 11.5 cm tip	0.5 mm
101-003-05	L= 11.5 cm tip	0.5 mm
101-008-05	L=11.5 cm tip	0.5 mm
101-006-10	L=12.1 cm tip	1.0 mm
101-016-10	L=16.0 cm tip	1.0 mm
101-018-10	L=18.0 cm tip	1.0 mm
101-020-10	L=20.0 cm tip	1.0 mm



## European Bipolar Forceps



Straight  
101-022-10



Bayonet  
101-036-10



Straight (Ang.)  
101-026-10



Bayonet  
101-038-10



Straight (Ang.)  
101-028-10



Straight (Ang.)  
101-030-10

Cat #	Discription	Size
101-022-10	L=22.0 cm tip	1.0 mm
101-026-10	L=16.0 cm tip	1.0 mm
101-028-10	L=18.0 cm tip	1.0 mm
101-030-10	L=20.0 cm tip	1.0 mm
101-036-10	L=16.0 cm tip	1.0 mm
101-038-10	L=18.0 cm tip	1.0 mm



## European Bipolar Forceps



Bayonet  
101-040-10



Bayonet (Ang.Up)  
101-046-10



Bayonet (Ang.Up)  
101-048-10



Bayonet (Ang.Up)  
101-050-10



Bayonet (Ang.Down)  
101-056-10



Bayonet (Ang.Up)  
101-058-10



Bayonet (Ang.Up)  
101-060-10

Cat #	Discription	Size
101-040-10	L=22.0 cm tip	1.0 mm
101-046-10	L=16.0 cm tip	1.0 mm
101-048-10	L=18.0 cm tip	1.0 mm
101-050-10	L=20.0 cm tip	1.0 mm
101-056-10	L=16.0 cm tip	1.0 mm
101-058-10	L=18.0 cm tip	1.0 mm
101-060-10	L=20.0 cm tip	1.0 mm





HIGH FREQUENCY INNOVATION

## European Irrigation Bipolar Forceps



Straight  
101-016-10 IR  
L=16.0 cm tip size: 1.0 mm



Straight  
101-020-10 IR  
L=20.0 cm tip size: 1.0 mm



Bayonet  
101-036-10 IR  
L=16.0 cm tip size: 1.0 mm



Bayonet  
101-040-10 IR  
L=20.0 cm tip size: 1.0 mm


Cat #	Discription	Size
101-016-10	L=16.0 cm tip	1.0 mm
101-020-10	L=20.0 cm tip	1.0 mm
101-036-10	L=16.0 cm tip	1.0 mm
101-040-10	L=20.0 cm tip	1.0 mm



## European Non-Stick Bipolar Forceps



McPherson (Str.)  
10-001-05 NS  
L= 9.0 cm tip size: 0.5 mm



McPherson (crv.)  
10-002-05 NS  
L= 11.5 cm tip size: 0.5 mm




Jeweler # 5  
10-003-05 NS  
L= 11.5 cm tip size: 0.5 mm



Adson  
10-006-10 NS  
L=12.1 cm tip size: 1.0 mm



Jeweler # 4  
10-007-05 NS  
L=11.5 cm tip size: 0.5 mm




Jeweler # 7  
10-008-05 NS  
L=11.5 cm tip size: 0.5 mm

Cat #	Discription	Size
10-001-05 NS	L= 9.0 cm tip	0.5 mm
10-002-05 NS	L= 11.5 cm tip	0.5 mm
10-003-05 NS	L= 11.5 cm tip	0.5 mm
10-006-10 NS	L=12.1 cm tip	1.0 mm
10-007-05 NS	L=11.5 cm tip	0.5 mm
10-008-05 NS	L=11.5 cm tip	0.5 mm



## European Non-Stick Bipolar Forceps



Coaptation  
10-009-05 NS

Coaptation Nadler  
10-010-05 NS

Iris (Str.)  
10-011-10 NS

Iris (crv.)  
10-012-10 NS  
size: 1.0 mm

Straight  
10-018-10 NS

Straight  
10-019-10 NS

Cat #	Description	Size
10-009-05 NS	L=10.8 cm tip	0.5 mm
10-010-05 NS	L=10.2 cm tip	0.5 mm
10-011-10 NS	L=10.8 cm tip	1.0 mm
10-012-10 NS	L=10.8 cm tip	1.0 mm
10-018-10 NS	L=18.0 cm tip	1.0 mm
10-019-10 NS	L=18.0 cm tip	1.0 mm





## European Non-Stick Bipolar Forceps



Straight  
10-020-10 NS



Straight  
10-022-10 NS  
size: 1.0 mm



Bayonet  
10-036-10 NS



Straight (Ang.)  
10-026-10 NS



Bayonet  
10-038-10 NS



Straight (Ang.)  
10-028-10 NS




Straight (Ang.)  
10-030-10 NS

Straight (Ang.)  
10-032-10 NS

Cat #	Discription	Size
10-020-10 NS	L=20.0 cm tip	1.0 mm
10-022-10 NS	L=22.0 cm tip	1.0 mm
10-026-10 NS	L=16.0 cm tip	1.0 mm
10-028-10 NS	L=18.0 cm tip	1.0 mm
10-030-10 NS	L=20.0 cm tip	1.0 mm
10-032-10 NS	L=22.0 cm tip	1.0 mm
10-036-10 NS	L=16.0 cm tip	1.0 mm
10-038-10 NS	L=18.0 cm tip	1.0 mm



## European Non-Stick Bipolar Forceps



Bayonet  
10-040-10 NS

Bayonet (Ang.Up)  
10-046-10 NS

Bayonet (Ang.Up)  
10-048-10 NS

Bayonet (Ang.Up)  
10-050-10 NS

Bayonet (Ang.Down)  
10-056-10 NS


Bayonet (Ang.Up)  
10-058-10 NS

Bayonet (Ang.Up)  
10-060-10 NS

Cat #	Discription	Size
10-040-10 NS	L=20.0 cm tip	1.0 mm
10-046-10 NS	L=16.0 cm tip	1.0 mm
10-048-10 NS	L=18.0 cm tip	1.0 mm
10-050-10 NS	L=20.0 cm tip	1.0 mm
10-056-10 NS	L=16.0 cm tip	1.0 mm
10-058-10 NS	L=18.0 cm tip	1.0 mm
10-060-10 NS	L=20.0 cm tip	1.0 mm



## USA 2 Pin Bipolar Forceps



Cat #	Discription	Size
08-001-05	L= 9.0 cm tip	0.5 mm
08-002-05	L=9.0 cm tip	0.5 mm
08-003-05	L=10.8 cm tip	0.5 mm
08-004-05	L=10.8 cm tip	0.5 mm
08-005-05	L=11.5 cm tip	0.5 mm
08-006-05	L=11.5 cm tip	0.5 mm
08-007-05	L=11.5 cm tip	0.5 mm
08-009-10	L=10.8 cm	1.0 mm
08-008-05	L=11.5 cm Tip	0.5 mm





## USA 2 Pin Bipolar Forceps



Cushing  
08-013-10



Cushing  
08-014-10



Cushing (Gursh)  
08-018-10



Iris (Crv.)  
08-010-10



Adson  
08-011-10



Semkin  
08-012-10

Cat #	Discription	Size
08-010-10	L=10.8 cm tip	1.0 mm
08-011-10	L=12.1 cm tip	1.0 mm
08-012-10	L=14.0 cm tip	1.0 mm
08-013-10	L=16.0 cm tip	1.0 mm
08-014-10	L=11.5 cm tip	1.0 mm
08--018-10	L=18.0 cm tip	1.0 mm



## USA 2 Pin Bipolar Forceps



**Semkin (Irrigation)**  
09-012IR



**Gerald Bayonet (Irrigation)**  
09-020IR



**Jansen Bayonet (Irrigation)**  
09-026IR



**Scollive Green-Wood (Irrigation)**  
09-031IR



**Hardy Bayonet (Irrigation)**  
09-037IR



## European Irrigation Bipolar Forceps

**Bayonet Forceps (Irrigation)**

11-036IR



**Bayonet Forceps (Irrigation)**

11-040IR



**Straight Forceps (Irrigation)**

11-020 IR



**Straight Forceps (Irrigation)**

11-014IR







## Disposable Bipolar Forceps



**Gerald Straight Single-Use**  
08-019S



**Jansen Bayonet (Irrigation)**  
08-028 IRS



**Jansen Bayonet Single-Use**  
08-030S



**Jeweler # 5 Single-Use**  
08-005S



**Adson Forceps Single-Use**  
08-011S



**Single-Use Jeweler # 5**

10-003S



**Single-Use Straight Forceps**

10-020S



**Single-Use Bayonet Forceps**

10-040S



**Single-Use Straight Forceps**

10-019S



**Single-Use Straight Forceps (ang.)**

10-030S



**Turnur Warwick Scissor (Straight)**  
05-001




**Turnur Warwick Scissor (Curved)**  
05-002




**Surgical Scissor (Straight)**  
05-011



**Surgical Scissor (Straight)**  
05-012




**Abbey Sub Needle**  
17.5cm 05-015



**Willson Hey Forceps (STR)**  
05-022



**Curved**  
05-024



**Artery Forceps (STR)**  
05-026



**Curved**  
05-26



**Richie Forceps (str)**  
05-038





**Adson**  
06-001



**Mcldone**  
06-002



**Gillie**  
06-001



**Disecting**  
06-001



**Gerald (1x2)**  
06-001



**Adson Bayonet**  
06-001



**Adson Bayonet (Ang. Up)**  
06-016



**Remington Hob**  
06-019



## Monopolar Cables



**Silicone Monopolar Cables 10" (3.0m)**  
 SI 19.230  
 Instrument Connector: Screw Type  
 Generator Connector: Standard Bovie 8.0mm Plug  
 Compatibility: Valleylab, Conmed, Bovie



**Silicone Monopolar Cables 10" (3.0m)**  
 SI 19.260  
 Instrument Connector: 4.0mm Male Fitting  
 Generator Connector: Standard 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed & All ES Generators



**Silicone Monopolar Endoscopic Cable 10" (3.0m)**  
 SI 19.240  
 Instrument Connector: 4.0mm Female Plug  
 Generator Connector: Standard Bovie 8.0mm Plug  
 Compatibility: Valleylab, Conmed, Bovie



**Silicone Monopolar Endoscopic Cable 10" (3.0m)**  
 SI 19.250  
 Instrument Connector: 4.0mm Female Plug  
 Generator Connector: Standard 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed & All ES Generators



**Silicone Diathermy Cable 10" (3.0m)**  
 SI 19.270  
 Instrument Connector: 4.8mm Male Fitting  
 Generator Connector: Standard Bovie 8.0mm Plug  
 Compatibility: Valleylab, Conmed, Bovie

**Silicone Diathermy Cable 10" (3.0m)**  
**SI 19.260**  
 Instrument Connector: 4.8mm Male Fitting  
 Generator Connector: Standard 4.0mm Banana Plug  
 Compatibility: Valleylab, Conmed & All ES Generators



SI 19.260



**Artery Forceps  
(Bi-Clamp) 18-004-25**

**Artery Forceps  
(Non Stick)  
18-006-25**

**Artery Forceps  
17-004-25**

**Artery Forceps  
17-001-25**



Silicone Bi-Clamp Connecting cord

13






## Bipolar Electrodes



**Straight**  
14-001




**Angled**  
14-002



**Angled**  
14-003



**Angled 45°**  
14-004




**Angled 90°**  
14-005



**Bi-Turbinate Probe (17.5cm)**  
14-028



## Electrosurgical Pencil



**Single-USE ESU Finger Switch Pencil**  
01-001



**Single-Use Foot Control Pencil**  
01-005



**Finger-Switch Pencil Single Use**  
01-004  
(3m Cable)  
(5m Cable)



**Reusable ESU Finger Switch Pencil**  
01-015



**Reusable ESU Foot Control Pencil**  
01-008

## Electrodes 2.4mm



Blade Electrode  
Disposable: 03-001



Ball Electrode (3.0 mm)  
Disposable: 03-007



Circle Loop Electrode 4.0 mm  
Reusable: 04-001



Ball Electrode (crv.) (3.0 mm)  
Disposable: 03-008



Blade Electrode (crv.)  
Disposable: 03-002



Circle Loop Electrode 4.0 mm  
Reusable: 04-001



Needle Electrode  
Disposable: 03-003



Ball Electrode (4.0 mm)  
Disposable: 03-009



Circle Loop Electrode 8.0 mm  
Reusable: 04-003



Needle Electrode (crv.)  
Disposable: 03-004



Ball Electrode (crv.) (4.0 mm)  
Disposable: 03-010



Diamond Electrode 8.0 mm  
Reusable: 04-002



Ball Electrode (2.0 mm)  
Disposable: 03-005



Ball Electrode (5.0 mm)  
Disposable: 04-011




Ball Electrode (crv.) (2.0 mm)  
Disposable: 03-006



Ball Electrode (crv.) (5.0 mm)  
Disposable: 03-012



Tip Cleaner 03-35



Holster 03-38



Blade Electrode (Shaft 13.0 cm)  
Disposable: 03-013




Blade Electrode (Shaft 13.0 cm)  
Disposable: 03-015



Blade Electrode (Shaft 13.0 cm)  
Disposable: 03-017



Blade Electrode (Shaft 13.0 cm)  
Disposable: 03-023



Coagulation Electrode  
Disposable: 03-025





BLADE ELECTRODE (Non-Stick)  
5.0 cm 03-51



BLADE ELECTRODE (Non-Stick)  
7.0 cm 03-52



BLADE ELECTRODE (Non-Stick)  
10.0 cm 03-53



BLADE ELECTRODE (Non-Stick)  
15.0 cm 03-54



Needle ELECTRODE (Non-Stick)  
7.0 cm 03-55



Needle ELECTRODE (Non-Stick)  
03-56 10.0 cm



Needle ELECTRODE (Non-Stick)  
03-57 15.0 cm



Ball ELECTRODE (Non-Stick)  
7.0 cm (3.0 mm)03-58



Blade ELECTRODE (Non-Stick)  
5.0 cm 03-61



Blade ELECTRODE (Non-Stick)  
7.0 cm 03-62



Blade ELECTRODE (Non-Stick)  
10.0 cm 03-63



Blade ELECTRODE (Non-Stick)  
15.0 cm 03-64



Needle ELECTRODE (Non-Stick)  
5.0 cm 03-65



Needle ELECTRODE (Non-Stick)  
10.0 cm 03-66



Needle ELECTRODE (Non-Stick)  
15.0 cm 03-67



Blade ELECTRODE (Non-Stick)  
15.0 cm (3.0 mm) 03-68



## Tungsten Needle Electrodes 2.4mm



Fine Needle ELECTRODE  
5.0 cm 03-81



Fine Needle ELECTRODE  
7.0 cm 03-82



Fine Needle ELECTRODE  
9.0 cm 03-83



Fine Needle ELECTRODE  
10.0 cm 03-84



Fine Needle ELECTRODE  
13.0 cm 03-85



Fine Needle ELECTRODE  
15.0 cm 03-86

"**Tungsten Needles** are widely used in a variety of surgical procedures like Plastic Surgery, ENT surgery, Head & Neck Surgery, Oral & Maxillofacial surgery. Tungsten has extremely high melting point (3400°C) which helps maintain a sharpness compared to stainless steel tips that dull rapidly"





## Tungsten Needle Electrodes 2.4mm



Fine Needle ELECTRODE (Crv.)  
5.0 cm 03-91



Fine Needle ELECTRODE (Crv.)  
7.0 cm 03-92



Fine Needle ELECTRODE (Crv.)  
9.0 cm 03-93



Fine Needle ELECTRODE (Crv.)  
10.0 cm 03-94



Fine Needle ELECTRODE (Crv.)  
13.0 cm 03-95



Fine Needle ELECTRODE (Crv.)  
15.0 cm 03-96

"**Tungsten Needles** are widely used in a variety of surgical procedures like Plastic Surgery, ENT surgery, Head & Neck Surgery, Oral & Maxillofacial surgery. Tungsten has extremely high melting point (3400°C) which helps maintain a sharpness compared to stainless steel tips that dull rapidly"



## Loop Electrodes 2.4mm



Size	1.0X0.4 cm	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm
Disposable	150-170-04	150-170-05	150-170-07	150-170-08	150-170-10
Reusable	152-170-04	152-170-05	152-170-07	152-170-08	152-170-10



Size	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm	1.0X1.5 cm
Disposable	150-175-05	150-175-07	150-175-08	150-175-10	150-175-15
Reusable	152-175-05	152-175-07	152-175-08	152-175-10	152-175-15



Size	1.5X0.5 cm	1.5X0.7 cm	1.5X0.8 cm	1.5X1.0 cm	1.5X1.2 cm	1.5X1.5 cm
Disposable	150-180-05	150-180-07	150-180-08	150-180-10	150-180-12	150-180-15
Reusable	152-180-05	152-180-07	152-180-08	152-180-10	152-180-12	152-180-15



Size	2.0X0.8 cm	2.0X1.0 cm	2.0X1.2 cm	2.0X1.3 cm	2.0X1.5 cm	2.0X2.0 cm
Disposable	150-190-08	150-190-10	150-190-12	150-190-13	150-190-15	150-190-20
Reusable	152-190-08	152-190-10	152-190-12	152-190-13	152-190-15	152-190-20



## Conization Electrodes



Conization Electrode  
Disposable 03-026




Conization Electrode  
Disposable 03-027



Conization Electrode  
Disposable 03-029



Conization Electrode  
Disposable 03-030



Conization Electrode  
Disposable 03-031



Extender Electrode 03-032  
7.0 cm



Extender Electrode 03-033  
13.0 cm



Extender Electrode 03-034  
18.0 cm





## Electrodes 4.0mm



**Blade Electrode**  
Reusable 20-001



**Ball Electrode**  
Reusable 20-007



**Blade Electrode (Crv)**  
Reusable 20-002



**Ball Electrode (Crv) 3mm**  
Reusable 20-008



**Needle Electrode**  
Reusable 20-003



**Ball Electrode (4mm)**  
Reusable 20-009



**Needle Electrode (Crv)**  
Reusable 20-004



**Ball Electrode (Crv) (4mm)**  
Reusable 20-010



**Ball Electrode (2mm)**  
Reusable 20-005



**Ball Electrode (5mm)**  
Reusable 20-011



**Ball Electrode (Crv) (2mm)**  
Reusable 20-006



**Ball Electrode (Crv) (5mm)**  
Reusable 20-012



**Blade Electrode (Shaft 13.0cm)**  
Reusable 20-018



**Needle Electrode (Shaft 13.0cm)**  
Reusable 20-019



**Ball Electrode 3mm (Shaft 13.0 cm)**  
Reusable 20-020



**Ball Electrode 4mm (Shaft 13.0cm)**  
Reusable 20-021



**Ball Electrode 5mm (Shaft 13.0cm)**  
Reusable 20-022



**Round Loop Electrode (7.0cm) 4.0mm**  
Reusable 20-013




**Round Loop Electrode (7.0cm) 6.0mm**  
Reusable 20-014




**Round Loop Electrode (7.0cm) 8.0mm**  
Reusable 20-015



**Round Loop Electrode (7.0cm) 10.0mm**  
Reusable 20-016



**Plate Electrode Small**  
Reusable 20-017



**Plate Electrode Large**  
Reusable 20-018



## Loop Electrodes 4.0mm



Size	1.0X0.4 cm	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm
Reusable	162-170-04	162-170-05	162-170-07	162-170-08	162-170-10



Size	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm	1.0X1.5 cm
Reusable	162-175-05	162-175-07	162-175-08	162-175-10	162-175-15




Size	1.5X0.5 cm	1.5X0.7 cm	1.5X0.8 cm	1.5X1.0 cm	1.5X1.2 cm	1.5X1.5 cm
Reusable	162-180-05	162-180-07	162-180-08	162-180-10	162-180-12	162-180-15




Size	2.0X0.8 cm	2.0X1.0 cm	2.0X1.2 cm	2.0X1.3 cm	2.0X1.5 cm	2.0X2.0 cm
Reusable	162-190-08	162-190-10	162-190-12	162-190-13	162-190-15	162-190-20



## Patient Return Plate



**Patient Return Plate Reusable**  
(5m Cable) 25-001



**Patient Return Plate Reusable**  
(5m Cable) 25-002



**Patient Return Cable**  
25-006-a (5m cable)  
25-006-b (3m Cable)

**Patient Return Cable**  
25-007-a (5m cable)  
25-007-b (3m Cable)



**S.S Patient Return Plate**  
(16x12cm)25-003






## Laparoscopic Electrodes


### WITHOUT SUCTION



15-001 L-Hook Electrode (36cm)




15-002 J-Hook Electrode (36cm)



15-003 Needle-Hook Electrode (36cm)



15-004 Spatula-Hook Electrode (36cm)



15-005 Ball-Hook Electrode (36cm)

### WITH SUCTION



15-006 L-Hook Electrode (Irrigation) (36cm)



15-007 J-Hook Electrode (Irrigation) (36cm)



15-008 Needle-Hook Electrode (Irrigation) (36cm)



15-009 Spatula-Hook Electrode (Irrigation) (36cm)



15-010 Ball-Hook Electrode (Irrigation) (36cm)

## Laparoscopic Electrodes






Dagger 17-100 (12cm)



Dagger 17-101 (18cm)




Hook 17-102 (12cm)



Hook 17-103 (18cm)




Hockey 17-104 (12 cm)



Hockey 17-105 (18 cm)



**Frazier Suction Tube**  
17-016  
11 French (5", 13.5cm)



**Frazier Suction Tube**  
17-016  
11 French (8", 20.5cm)



**Frazier Suction Tube**  
17-016  
11 French (16", 40.6cm)



## Electrosurgical Instruments for Gynecology



**Graves Vaginal Speculum**  
with smoke evacuation tube  
Small: 21-001  
Medium: 21-002  
Large: 21-003



**VU-Max Speculum**  
with smoke evacuation tube  
Small: 21-011  
Medium: 21-012  
Large: 21-013



**Pederson Vaginal Speculum**  
with smoke evacuation tube  
Small: 21-041  
Medium: 21-042  
Large: 21-043



**VU-More (Graves) Speculum**  
with smoke evacuation tube  
Small: 21-031  
Medium: 21-032  
Large: 21-033



**Collin Speculum**  
with smoke evacuation tube  
Small: 21-051  
Medium: 21-052  
Large: 21-053



**CUSCO Speculum**  
with smoke evacuation tube  
Small: 21-021  
Medium: 21-022  
Large: 21-023





## Electrosurgical Instruments for Gynecology



**Weisman Speculum  
with Suction (Left Open)**  
Small: 21-061  
Medium: 21-062  
Large: 21-063



**Weisman Speculum  
with Suction (Right Open)**  
Small: 21-071  
Medium: 21-072  
Large: 21-073

**Dressing Forceps**  
23-001 (21.0 cm)  
23-001 (25.0 cm)



**Tissue Forceps**  
23-003 (21.0 cm)  
23-004 (25.0 cm)



**Iris Hook**  
23-005 (25.0 cm)



**Emmet Hook**  
23-006 (25.0 cm)



**Straight Hook**  
23-007 (25.0 cm)



**Two-Prong Hook**  
23-008 (25.0 cm)



**Three-Prong Hook**  
23-009 (25.0 cm)






## Electrosurgical Instruments for Gynecology



**Lateral Vaginal  
Retractor**  
22-001



**Kogan Endocervical  
Speculum**  
22-003



**Lateral Vaginal  
Retractor (wide-Open)**  
22-006



**Rigby Retractor**  
22-007



**Rigby Retractor**  
22-005








## Electrosurgical Instruments for Gynecology



**Schroeder Tanaculum**  
22-008 (24.0 cm)




**Schroeder Tanaculum**  
24-009 (24.0 cm)



**Single-Tooth Tanaculum**  
22-010 (24.0 cm)



**Champion Forceps**  
22-012 (22.0 cm)



**Ring Forceps**  
22-011 (22.0 cm)





## Electrosurgical Instruments for Gynecology



**Kogan Endocervical Speculum**  
25-003



**Lateral Vaginal Retractor**  
25-001




**Graves Speculum**  
with smoking evacuation tube  
Small: 21-001  
Medium: 21-002  
Large: 21-003



**Pederson Graves Speculum**  
with smoking evacuation tube  
Small: 21-001  
Medium: 21-002  
Large: 21-003



**Collin Speculum**  
with smoking evacuation tube  
Small: 24-051  
Medium: 24-052  
Large: 24-053






**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.100**

Instrument Connector: USA 2-Pin Plug  
Generator Connector: 4.0mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olpymus and Mbst Bipolar Generator



**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.120**

Instrument Connector: European Flat Plug  
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Wolf, Berchold



**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.110**

Instrument Connector: European Flat Plug  
Generator Connector: 4.0mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olpymus & Mbst Bipolar Generator



**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.130**

Instrument Connector: USA2-Pin Plug  
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Wolf, Berchold



**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.150**

Instrument Connector: USA2-Pin Plug  
Generator Connector: Erbe Type


Compatibility: Eibe, Siemens, Downs



**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.170**

Instrument Connector: USA2-Pin Plug  
Generator Connector: Pin Distance 29.0mm

Compatibility: Valleylab, EMS, Dolley



**Silicone Bipolar Cables 10'' (3.0m)**  
**SI 19.160**

Instrument Connector: European Flat Plug  
Generator Connector: Pin Distance 29.0mm

Compatibility: Valleylab, EMS, Dolley



# Surgimpex & Co.

Street Sadat School Shahab Pura  
Sialkot-51310 PAKISTAN.

Cell: +92 331 6163388 / Tel: +92 52 3563388  
Email: [surgimpexco@gmail.com](mailto:surgimpexco@gmail.com) / [info@surgimpex.com](mailto:info@surgimpex.com)  
URL: [www.surgimpex.com](http://www.surgimpex.com)