

# Surgimpex & Co.

Laryngoscope, Otoscope, Ophthalmoscope



## Diagnostic Instruments



ایران صنعت و تجارت سیالکوٹ

**History** To introduce the firm briefly , M/s. Surgimpex & Co was established in 1978 Mr. Abdul Hafiz Mughal , who was well known as a specialized skilled worker and maker of Surgical Dental Veterinary Instruments.

### **Our Quality**

1 Fair dealing, Straight Forwardness, Honesty, Sincerity, Truth

2 The Material used in the production against order is exactly as per approved samples.

3 Normally, in Cutting Instruments, we use (420 Type) Japanese Stainless Steel.

4 For non -Cutting instruments like Forceps, Tweezers etc.; we use (410 Type) Japanese Stainless Steel.

5 Moreover, we can also the material as per our customer requirements.

### **Our Price**

1 To provide “Excellent Service and Economically Price” is the motto of our business

2 All instruments are made of Stainless Steel otherwise specified.

3 “All types of Mirror-Finish” “Satin/Dull-Finish” “Sand Blast-Finish” instruments are available in stock but choice of polishes will be supplied on demand without any extra charges

Therefore, You are requested to kindly provide us a list of your desire items against which we shall appreciate to quote / offer our most competitive prices for favour of your kind approval.

Looking forward to your early kind encouraging reply.

### **Our Services**

#### **Manufacturing Steps**

#### **Production Capabilities**

Surgimpex & Co is an employee owned company of craftsmen and has got a good reputation for quality and reliability. Our significant manufacturing steps are being described here-under briefly

1.MATERIAL USED: Always a best quality stainless steel material is purchased from steel manufacturers. Two kind of series i.e. 300, 410 and 420 of stainless steel is used for production of Surgical Instruments.

mechanization to production run



# Reusable Laryngoscope



LED/ Halogen Bulb

## MACINTOSH

Integrated flexible Laryngoscope (Fiber Optic) Cold light Guide.

MACINTOSH	Fig	Working Length
SI - 1001-01	1	72 mm
SI - 1001-02	2	88 mm
SI - 1001-03	3	112 mm
SI - 1001-04	4	134 mm

## FEATURES

303/304 Surgical Stainless Steel  
Smooth edges. Rounded tip and high quality surface finish large fiber blades offers more light no black mark fully autoclavable and easy to clean stain finish.

## MACINTOSH

Flexible interchangeable blade (fiber optic) cold light guide.

MACINTOSH	Fig	Working Length
SI - 1003-02	2	80 mm
SI - 1003-03	3	92 mm
SI - 1003-04	4	116 mm
SI - 1003-05	5	135 mm

## FEATURES:

303/304 Surgical stainless steel  
Smooth edges. Rounded tip and high quality surface finish removable light guides design allow easier cleaning and extends blade life provides with 4mm fiber optic bundles or light pass gives true tissue color.



LED/ Halogen Bulb



## MACINTOSH

Flexible fix tube fiber optic cold light guide.

MACINTOSH	Fig	Working Length
SI - 1005-01	1	70 mm
SI - 1005-02	2	88 mm
SI - 1005-03	3	112 mm
SI - 1005-04	4	134 mm

## FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. Rounded tip and high quality surface finish provides with 3mm Fiber Optic  
Offers no black marks fully autoclavable and easier to clean stain finish.



LED/ Halogen Bulb



LED/ Halogen Bulb

## MACINTOSH

Flexible fix tube fiber optic cold light guide.

MACINTOSH	Fig	Working Length
SI - 1007-02	2	92 mm
SI - 1007-03	3	113 mm
SI - 1007-04	4	135 mm
SI - 1007-05	5	155 mm

## FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. Rounded tip and high quality surface finish. Three spring ball bearings gives secure fit and firm support supplied with LED or halogen bulb satin finish.



## Integrated Laryngoscope (Fiber Optic) Cold Light Guide

MACINTOSH	Fig	Working Length
SI - 1009-01	1	72 mm
SI - 1009-02	2	90 mm
SI - 1009-03	3	112 mm
SI - 1009-04	4	135 mm



### FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. Rounded tip and high quality surface finish large Fiber bundles offers more light no black marks fully Autocleavable and easier to clean satin finish.




LED/ Halogen Bulb



LED/ Halogen Bulb

### FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. high quality surface finish provides with 3mm Fiber Optic offers no Black Marks fully Autoclavable and easier to clean satin finish.



## Miller Fix Tube HBT (Fiber Optic) Cold Light Guide

MILLER	Fig	Working Length
SI - 1011-00	00	45 mm
SI - 1011-0	0	55 mm
SI - 1011-01	1	80 mm
SI - 1011-02	2	130 mm
SI - 1011-03	3	170 mm
SI - 1011-04	4	180 mm






LED/ Halogen Bulb

### Macintosh Interchangeable (Fiber Optic) Cold Light Guide

MACINTOSH	Fig	Working Length
SI - 1013-01	1	80 mm
SI - 1013-02	2	92 mm
SI - 1013-03	3	116 mm
SI - 1013-04	4	135 mm

### FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. Rounded tip and high quality surface finish. Removable Light guides design allow easier cleaning and extends blade life provides with 4mm Fiber Optic bundles for light pass gives true tissue color satin finish.



### Miller Fix tube HBT (Fiber Optic) Cold Light Guide.

MILLER	Fig	Working Length
SI - 1015-00	00	45 mm
SI - 1015-0	0	55 mm
SI - 1015-01	1	80 mm
SI - 1015-02	2	130 mm
SI - 1015-03	3	170 mm
SI - 1015-04	4	180 mm

### FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. High quality surface finish. Provides with 4mm Fiber Optic Bundles for light pass gives true tissue color fully autoclavable and easier to clean satin finish



LED/ Halogen Bulb



## Macintosh Fix Tube (Fiber Optic) Cold Light Guide


MACINTOSH	Fig	Working Length
SI - 1017-01	1	70 mm
SI - 1017-02	2	88 mm
SI - 1017-03	3	112 mm
SI - 1017-04	4	134 mm

### FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. Rounded tip and High quality surface finish. Provides with 3mm Fiber Optic Bundles or light pass gives true tissue color fully autoclavable and easier to clean satin finish.



LED/ Halogen Bulb



## Miller Fix Tube (Fiber Optic) Cold Light Guide

MILLER	Fig	Working Length
SI - 1019-00	00	45 mm
SI - 1019-0	0	57 mm
SI - 1019-01	1	80 mm
SI - 1019-02	2	130 mm
SI - 1019-03	3	170 mm
SI - 1019-04	4	180 mm

### FEATURES

303/304 Surgical Stainless Steel  
Smooth Edges. High quality surface finish. Provides with 3mm Fiber Optic Bundles for light pass gives true tissue color fully autoclavable and easier to clean satin finish.



LED/ Halogen Bulb




### FEATURES

303/304 Surgical Stainless Steel  
 Smooth Edges. Rounded tip and high quality surface finish. Three Spring ball bearing gives secure fit and firm support supplied with LED or Halogen Bulb  
 Satin finish.

### Macintosh Conventional

MACINTOSH	Fig	Working Length
SI - 1021-00	00	60 mm
SI - 1021-0	0	68 mm
SI - 1021-01	1	92 mm
SI - 1021-02	2	113 mm
SI - 1021-03	3	135 mm
SI - 1021-04	4	155 mm

LED/ Halogen Bulb



### Miller Conventional

MILLER	Fig	Working Length
SI - 1023-00	00	60 mm
SI - 1023-0	0	70 mm
SI - 1023-01	1	92 mm
SI - 1023-02	2	113 mm
SI - 1023-03	3	135 mm
SI - 1023-04	4	155 mm


### FEATURES

303/304 Surgical Stainless Steel  
 Smooth Edges. High quality surface finish.  
 Three spring ball bearing gives secure fit and firm support supplied with LED or Halogen Bulb  
 Satin finish.



LED/ Halogen Bulb







### Macintosh Conventional co

MACINTOSH	Fig	Working Length
SI - 1025-01	1	80 mm
SI - 1025-02	2	92 mm
SI - 1025-03	3	116 mm
SI - 1025-04	4	135 mm

### FEATURES

303/304 Surgical Stainless Steel  
 Smooth Edges. Rounded tip and High quality surface finish. Three spring ball bearings gives secure fit and firm support supplied with LED 2.5v 5000 Lux satin finish



### Macintosh Conventional co

MILLER	Fig	Working Length
SI - 1027-00	00	45 mm
SI - 1027-0	0	57 mm
SI - 1027-01	1	80 mm
SI - 1027-02	2	130 mm
SI - 1027-03	3	170 mm
SI - 1027-04	4	180 mm

### FEATURES

303/304 Surgical Stainless Steel  
 Smooth Edges. High quality surface finish.  
 Three spring ball bearing gives secure fit and firm support supplied with LED 2.5v 5000 Lux satin finish



**Macintosh F/O Fix Tube (GS 100)**  
Cold Light Guide

MACINTOSH	Fig	Working Length
SI - 1029-01	1	80 mm
SI - 1029-02	2	92 mm
SI - 1029-03	3	116 mm
SI - 1029-04	4	135 mm



**FEATURES**

303/304 Surgical Stainless Steel  
Smooth Edges. Rounded tip and high quality surface finish. Provides with 3mm Fiber Optic offer no black marks fully autoclavable and easier to clean finish



LED/ Halogen Bulb

## Single USE Laryngoscope

**Macintosh**

Flexible conventional blade

MACINTOSH	Fig	Working Length
SI - 1031-01	1	76 mm
SI - 1031-02	2	88 mm
SI - 1031-03	3	112 mm
SI - 1031-04	4	130 mm

**FEATURES**

Premium quality Surgical Stainless Steel  
Spatula provides with unbreakable plastic base Unique heal design prevents Blade from touching the handle clearly  
Marked and packed single patient use  
smooth edges, rounded tip and high quality surface finish low cost than reusable one.



LED/ Halogen Bulb

# Single USE Laryngoscope



## MACINTOSH

Flexible Blade Fiber Optic

MACINTOSH	Fig	Working Length
SI - 1033-01	1	76 mm
SI - 1033-02	2	88 mm
SI - 1033-03	3	112 mm
SI - 1033-04	4	130 mm

## FEATURES

Premium quality surgical stainless steel spatula provides with unbreakable plastic base unique heal design p vents blade form touching the handle clearly marked and packed single patie t use Fiber Optic tube black jac et end lit fiber smooth edges, rounded tip and high quality surface finish low cost reusable one.

## Flexible Conventional Blade

MACINTOSH	Fig	Working Length
SI - 1035-01	1	76 mm
SI - 1035-02	2	88 mm
SI - 1035-03	3	112 mm
SI - 1035-04	4	130 mm

## FEATURES

Premium quality surgical stainless steel spatula provides with unbreakable plastic base unique heal design prevents blade from touching the handle clearly marked and packed single patie t use smooth edges, rounded tip and high quality surface finish low cost reusable one



LED Bulb





### Macintosh F/O Fix Single Use

MACINTOSH	Fig	Working Length
SI - 1037-01	1	82 mm
SI - 1037-02	2	92 mm
SI - 1037-03	3	116 mm
SI - 1037-04	4	135 mm

### Features

Tip to heal 303/304 Surgical Stainless steel  
Smooth edges, Rounded tip and high quality surface finish.

Three spring ball bearings gives secure fit and firm support Fiber Optic tube back jacket and Lit fiber, satin finish.



LED/Halogen Bulb



### Macintosh Eco Single Use

MACINTOSH	Fig	Working Length
SI - 1039-01	1	82 mm
SI - 1039-02	2	92 mm
SI - 1039-03	3	116 mm
SI - 1039-04	4	135 mm

### Features

Tip to heal 303/304 Surgical stainless steel  
smooth edges, rounded tip and high quality surface finish.

Three spring ball bearings gives secure fit and firm support plastic sleeve 2.5v LED 5000 Lux Brightness satin finish.





## Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
SI - 1041-01	1	79 mm
SI - 1041-02	2	88 mm
SI - 1041-03	3	112 mm
SI - 1041-04	4	130 mm

### FEATURES

Tip to Spatula 303/304 Surgical Stainless Steel  
 Heel silver smooth edges, Rounded tip and  
 High quality surface finish.  
 Three spring ball bearings gives secure fit and  
 firm support fiber optic tube black jet and  
 LIT Fiber, satin finish



## Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
SI - 1043-01	1	76 mm
SI - 1043-02	2	88 mm
SI - 1043-03	3	112 mm
SI - 1043-04	4	130 mm

### Features

Tip to Spatula 303/304 Surgical Stainless Steel  
 Heel silver smooth edges, Rounded tip and  
 High quality surface finish.  
 Three spring ball bearings gives secure fit and  
 firm support fiber optic Acrylic and satin finish






**Macintosh Single Use (Fiber Optic)**

MACINTOSH	Fig	Working Length
SI - 1045-01	1	76 mm
SI - 1045-02	2	88 mm
SI - 1045-03	3	112 mm
SI - 1045-04	4	130 mm

**Features**

Tip to Spatula 303/304 Surgical Stainless steel  
 heel plastic (ABS) smooth edges, rounded tip and  
 high quality surface finish.

Two spring ball bearings gives secure fit and  
 firm support fiber optic tube black jet and  
 LIT Fiber, satin finish



**Miller Single Use (Fiber Optic)**


MILLER	Fig	Working Length
SI - 1047-00	00	45 mm
SI - 1047-0	0	55 mm
SI - 1047-01	1	80 mm
SI - 1047-02	2	130 mm
SI - 1047-03	3	170 mm
SI - 1047-04	4	180 mm

**Features**

Spatula 303/304 Surgical Stainless steel  
 Heel plastic (ABS) smooth edges,  
 High quality surface finish  
 Two spring ball bearings gives secure fit and  
 firm support fiber optic tube black jet end  
 Lit fiber , satin finish








## Macintosh Single Use (Fiber Optic)

MACINTOSH	Fig	Working Length
SI - 1049-01	1	76 mm
SI - 1049-02	2	88 mm
SI - 1049-03	3	112 mm
SI - 1049-04	4	130 mm

### FEATURES

Tip to Spatula 303/304 Surgical Stainless Steel  
 Heel plastic (ABS) smooth edges, rounded tip and high quality surface finish.  
 Two Spring ball bearings gives secure fit and firm support fiber acrylic rod satin finish




## Macintosh Conventional Single Use

MACINTOSH	Fig	Working Length
SI - 1051-01	1	83 mm
SI - 1051-02	2	95 mm
SI - 1051-03	3	120 mm
SI - 1051-04	4	138 mm

### Features

Tip to Spatula 303/304 Surgical Stainless Steel  
 Unique Heel design smooth edges, rounded tip and high quality surface finish.  
 Two spring ball bearings gives secure fit and firm support 2.5 single and double LED Design satin finish





### Miller Conventional Single Use

MILLER	Fig	Working Length
SI - 1053-00	00	45 mm
SI - 1053-0	0	55 mm
SI - 1053-01	1	82 mm
SI - 1053-02	2	130 mm
SI - 1053-03	3	135 mm
SI - 1053-04	4	185 mm

### Features

Tip to Spatula 303/304 Surgical Stainless Steel  
 Unique heel design smooth edges,  
 High Quality surface finish  
 Two spring ball bearings gives secure fit and firm support 2.5 single and double LED Design satin finish



### Macintosh Conventional Single Use

MACINTOSH	Fig	Working Length
SI - 1055-01	1	80 mm
SI - 1055-02	2	92 mm
SI - 1055-03	3	116 mm
SI - 1055-04	4	133 mm

### Features

Tip to Spatula 303/304 Surgical Stainless Steel  
 Unique Heel Design Smooth Edges, Rounded tip and  
 High Quality Surface Finish  
 Two Spring ball hearings gives secure fit and firm support 2.5 Single LED Design 5000 Lux Brightness  
 Satin Finish







## Conventional Sandad



Medium Small Stubby  
SI - 1057-01

## Fiber Standard



Medium Small Stubby  
SI - 1057-02

## Single Use Handles



Conventional Conventional Fiber Optic  
SI - 1057-03

## Hard Box & Soft pouches



SI - 1059-01



SI - 1059-02




SI - 1059-03



SI - 1059-04





SI - 1060-01



SI - 1061-01



SI - 1061-02



SI - 1062-01



SI - 1063-01



SI - 1063-02

# METAL OTOSCOPES



### FEATURES

Equipped with an easy to use Bayonet Connecting stem with high quality Professional metal head and chrome plated knurled battery handle making the surface a secure grip in the hand with extra bright light of 2.5v LED bulb illumination provided through Fiber Optic stem. Three times high quality magnifying lens.  
**Requires:** 20xC Batteries (Not included with set)  
**Packing:** In Black soft zipper pouch with 10 specula having Dia 2.5, & 4.0mm

**Fiber Optic  
Conventional Ooscope**  
SI-1065-01



### FEATURES

Equipped with an easy to use Bayonet Connecting stem with high quality Professional metal head and chrome plated knurled battery handle making the surface a secure grip in the hand with extra bright light of 2.5v LED bulb illumination provided through Fiber Optic stem. Three times high quality magnifying lens.  
**Requires:** 20xC Batteries (Not included with set)  
**Packing:** In Black soft zipper pouch with 10 specula having Dia 2.5, & 4.0mm

**Fiber Optic  
Metal Otoscope**  
SI-1066-01



### FEATURES

Standard Metal Otoscope with bayonet connecting head fitting ide or pneumaumatic tympanic membrane test and very day ENT use. 2.5v LED white light. Available in satin and polish finish.  
**Requires:** 2xC Batteries (Not Included with set)  
**Packing:** in black soft zipper pouch with 10 specula Having Dia 2.5mm & 4.0mm

**Conventional  
Metal Otoscope**  
SI-1067-01



### FEATURES

High Quality Professional Parker Otoscope with extra bright light 2.5v LED Illumination Three times high quality magnifying lens. Available in satin and polish finish.  
**Requires:** 2xC Batteries (Not included with set)  
**Packing:** in black soft zipper pouch with 10 specula having Dia 2.5mm & 4.0mm

**Conventional  
Parker Metal Otoscope**  
SI-1068-01





**FEATURES**

High Quality Professional Parker Otoscope with screw fitting illumination t. 20.5v LED Metal Handle & Plastic head CP . Rheostat dimmer

**Requires:** 2xC Bateriaes (Not included with set)  
**Packing:** In black soft zipper pouches with 10 Specula having Dia 2.5mm & 4.0mm

Conventional Parker Metal Otoscope  
SI-1070-01



**FEATURES**

2.5v Standard illumination LED Bulb. Chrome Plated Handle and 4x Magnifying Glasses. Reusable Speculum Size 0, 2.5, 3.5 and 4mm include. Connecting or pneumatic otoscope battery handle with rheostat for adjustment of illumination.

**Supply:** German bulb technology for super view with 3 reusable ear specula.  
**Requires:** 2xC Bateriaes (Not included with set)  
**Packing:** In black soft zipper pouches Available in satin and polish finish.

Conventional Metal Otoscope  
SI-1071-01



**FEATURES**

Open System with rotatable lens and optional speculum for convenient procedures 2.5v Otoscope.

**Requires:** 2xC Bateriaes (Not included with set)  
**Packing:** In black soft zipper pouches with 3 VET reusable specula having Dia 4.0, 5.0, 7.0mm

Conventional Veterinary Metal Otoscope  
SI-1072-01



**FEATURES**

Having special focal distance with two fold magnification. yonet closure. 2.5v Otoscope.

**Requires:** 2xC Bateriaes (Not included with set)  
**Packing:** In black soft zipper pouches with 3 VET reusable specula having Dia 4.0, 5.0, 7.0mm

Conventional Veterinary Metal Otoscope  
SI-1073-01



# MINI OTOSCOPES



### FEATURES

High Quality Pocket Sized mini Fiber Optic Otoscope with extra bright light 2.5v with 3x 4mm bundle with enhanced LED technology.  
**Requires:** 2x AA Bateriaes (Not included with set)  
**Packing:** In Gray or black pouch with 10 disposable specula Dia2.5mm & 40.0mm

**Available Colors:** ■ ■ ■ ■ ■ ■ ■ ■

**Fiber Optic  
Pocket Style**  
SI-1075-01



**Fiber Optic  
Pocket Style**  
SI-1076-01



### FEATURES

High Quality Pocket Sized mini Fiber Optic Otoscope with extra bright light 2.5v with 3x Magnification lens. Dal fi fivivver optic illu mination 4mm bundle with enhanced LED technology.  
**Requires:** 2x AA Bateriaes (Not included with set)  
**Packing:** In Gray or black pouch with 10 disposable specula Dia2.5mm & 40.0mm

**Available Colors:** ■ ■ ■ ■ ■ ■ ■ ■

**Conventional  
Pocket Style**  
SI-1077-01



**Conventional  
Pocket Style**  
SI-1078-01



### FEATURES

Top Quality Ophthalmoscope with fixation aperture in pocket format. Quality optic in miniature, Ideal for rounds and house calls. Attachment clip with integrated on/off switch. Secure, switches automatically when returned to the pocket. High quality handle. White upper section finished plastic shockproof, sturdy, non-slip. Enhanced LED Light Technology concentrated white light for perfect illumination and a brilliant image. 4 different apertures, A full-function instrument including fixation, Range of lenses, 20D to +20D, 2-piece handle and head, Easy to maintain, flexible

**Replaceable Batteries Size AA**  
**Supply:** Bristle and Hard Plastic Case Packing

**Ophthalmoscope**  
SI-1080-01



The Ophthalmoscope is a tool used in medicine to examine the interior of the eye, studying the inner parts of the eye, such as the optic disk, retinal blood vessels, retinal choroid and macula, can help diagnose diseases of the eyes.

**Ophthalmoscope**  
SI-1081-01



SI-1082-01  
**FEATURES**

GM Ophthalmoscope set pin contacting fitting ophthalmoscope head pin contact fitting lens selection wheel permits choice of correction lenses from 0 to -25 and 0 to +40 diameters single large circle aperture, rheostat for adjustment of illumination 2 battery handle choice

**Require:** 2x C Size batteries (Not Included with sets)  
**Supply:** Hard Case



SI-1083-01  
**FEATURES**

Parker mini ophthalmoscope 2.5v LED Bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses -25 to +40 for excellent resolution (red for - Green for +) Rubber eye brow rest to prevent scratching for eye-glasses illumination lens dial or clear identification markings or adaption settings.

**Require:** 2 x AA Size batteries (Not Included with sets)  
**Supply:** Soft pouch or plastic hard case.



SI-1084-01  
**FEATURES**

Parker mini ophthalmoscope 2.5v LED Bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses -25 to +40 for excellent resolution (red for - Green for +) Rubber eye brow rest to prevent scratching for eye-glasses illumination lens dial or clear identification markings or adaption settings.

**Require:** 2 x C Size batteries (Not Included with sets)  
**Supply:** Soft pouch or plastic hard case.

# Surgimpex & Co.

Laryngoscope

Otoscope

Ophthalmoscope



E-MAIL

Info@surgimpex.com



PHONE

+923316163388



WEBSITE

www.surgimpex.com



ADDRESS

Defence Road Mustafa Town

S I E Sialkot Pakistan

